

Hanzl® HanoPlus® 1L

Skin care cream with natural oils, betaine, bisabolol and beeswax



Description:	HANZL® HANOPLUS® is a white, slightly greasy and well absorbed cream with natural oils, betaine, bisabolol and beeswax. Silicone-free.
Composition {INCI}:	AQUA (WATER), ETHYLHEXYL STEARATE, GLYCERYL STEARATE SE, STEARYL ALCOHOL, PROPYLENE GLYCOL, GLYCERIN, THEOBROMA CACAO (COCOA) SEED BUTTER, CAPRYLIC/CAPRIC TRIGLYCERIDE, CERA ALBA (BEESWAX), BUTYROSPERMUM PARKII (SHEA) BUTTER, BETAINE, BISABOLOL, PERSEA GRATISSIMA (AVOCADO) OIL, PRUNUS AMYGDALUS DULCIS (SWEET ALMOND) OIL, CARBOMER, SODIUM CETEARYL SULFATE, PHENYLPROPANOL, PARFUM (FRAGRANCE), TOCOPHEROL, PHENOXYETHANOL, PROPANEDIOL, CAPRYLYL GLYCOL, P-ANISIC ACID, SODIUM HYDROXIDE
Properties:	HANZL® HANOPLUS® is a skin care cream for hands and face. Specifically selected natural oils and lipids such as almond oil, avocado oil and cocoa butter help keeping the skin soft and supple, without leaving a "greasy" feeling. The 3B skin care complex of betaine, bisabolol and beeswax ensures effective skin care. Regular use helps to keep the skin healthy and strengthens its natural defences.
Fields of application:	HANZL® HANOPLUS® is intended for care of skins under both professional and environmental stress. HANZL® HANOPLUS® is mainly used to regenerate the skin after work or during the breaks and is especially recommended for dry, rough and chapped skin. HANZL® HANOPLUS® is suited for use on hands and face.
Use:	Apply HANZL® HANOPLUS® to clean, dry skin after work as well as in breaks. Rub in thoroughly, paying particular attention to between the fingers and around the nails.
Regulatory framework:	HANZL® HANOPLUS® is subject to the EU Cosmetics Regulations and corresponds to it in its current version. It does not have to be labelled according to the regulations on hazardous substances.
Storage and shelf life:	HANZL® HANOPLUS® can be kept for 30 months at room temperature (from date of manufacture). The date of expiry is denoted on the packaging.
Packaging size:	1000 ml pouch bottle Dispensing-system: DUROL!NE® 1L dispenser.
Disposal:	The cardboard boxes, tubes and plastic packaging are recyclable.